Niti-S & ComVi Biliary Stent

Authorized representative in Europe
Représentant autorisé en Europe
Bevollmächtigter in Europa
Rappresentante autorizzato per l'Europa
Representante autorizado en Europa
Avrupa yetkili temsilcisi
Autoriseret repræsentant i Europa
유럽 대리인

Consult instructions for use
Consulter les instructions d'utilisation
Gebrauchsanleitung beachten
Consultare le istruzioni per l'uso
Consulte las instrucciones de uso
Kullanım talimatlarına bakın
Konsultér brugsanvisning

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Catalogue No.
N° de référence
Katalog-Nr.
N. di catalogo
N° de catálogo
Katalog No.
Referenzenummer
형명

Attention, consult instructions for use
Attention, consulter les instructions d'utilisation
Achtung: Gebrauchsanleitung beachten
Attenzione, consultare le istruzioni per l'uso
Atención, consulte las instrucciones de uso
Dikkat! Kullanım talimatlarına bakın.
Vær opmærksom, konsultér brugsanvisning
주의,사용상 지시 참조

Serial No.
N° de série
Serien-Nr.
N. Seriale
N° de serie
Seri No.
Serienummer
일련번호

Date of Manufacture
Date de fabrication
Datum der Herstellung
Data di produzione
Fecha de fabricación
Üretim Tarihi
Fabrikationsdato
제조일자

Do not use if package is damaged
Ne pas utiliser si l'emballage est abîmé
Nicht verwenden, wenn die Verpackung beschädigt ist!
Non utilizzare la confezione se danneggiata
No usar si el paquete está dañado.
Ambalaj hasar görmüşse tekrar kullanmayın
Må ikke bruges, hvis emballagen er beskadiget
포장이 파손된 경우 사용금지

Temperature limitation
Limites de température
Temperaturvorgaben
Intervallo termico
Límites de temperatura
Isi sınırlaması
Temperaturgrænser
보관 온도

Manufacturer
Fabricant
Hersteller
Produttore
Fabricante
Üretici
Fabrikant
제조사

Sterilized using ethylene oxide
Stérilisé à l'aide d'oxyde d'éthylène
Mit Ethylenoxid sterilisiert
Sterilizzato con ossido di etilene
Esterilizado por óxido de etileno
Etilen oksit kullanılarak sterilize edilmiştir
Steriliseret med ethylenoxid
E.O Gas 멸균

Use by(Expiration Date)
À utiliser avant (date d'expiration)
Verwendbar bis (Verfallsdatum)
Utilizzare entro (data di scadenza)
Utilizar antes de (fecha de vencimiento)
Son Kullanma Tarihi
Bruges før (udløbsdato)
유효기간

Do not reuse
Ne pas réutiliser
Nicht wiederverwenden
Non riutilizzare
No reutilizar
Yeniden kullanmayın
Må ikke genbruges
재사용 금지

Do not resterilize
Ne pas stériliser à nouveau
Nur einmal sterilisieren!
Non risterilizzare
No reesterilizar
Yeniden sterilize etmeyin
Må ikke reesteriliseres
재멸균 금지

MR Conditional
RM conditionnelle
Bedingte magnetische Resonanzsicherheit
RM compatibile
RM condicional
Koşullu MR
MR-betinget
MR 조건부 입증

TaeWoong
MEDICAL

MRI Information

MR Conditional

Niti-S & comvi stent was determined to be MR-conditional. Non-clinical testing demonstrated that the Niti-S & comvi stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, Niti-S & comvi stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position Niti-S & comvi stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm ²	676-mm ²	7,984-mm ²	977-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

User's Manual

1. Description

The Niti-S & ComVi Biliary Stent consists of the implantable metallic stent and introducer system.

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Niti-S Biliary Uncovered Stent
Niti-S Biliary Covered Stent
ComVi Biliary Stent

Figure 1. Stent Model

The Stent is loaded in introducer system and upon deployment the stent imparts an outward radial force on the luminal surface of the bile duct to establish patency.

Niti-S Full Covered Biliary Stents used in benign stricture can be removed; (see Warnings).

Full Covered Biliary Stents can be repositioned after deployment (see Warnings).


Figure 2. Introducer System (Percutaneous & Endoscopic)

- The percutaneous introducer system has a usable length of 50cm
- The endoscopic introducer system has a usable length of 180cm

Percutaneous Type is recommended

- When approached percutaneously

Endoscopic Type is recommended

- When approached endoscopically

2. Principle of Operation

The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher towards the hub.

Retraction of the outer sheath releases the stent.

3. Indication for Use

The Niti-S & ComVi Biliary Stent is intended for maintaining biliary luminal patency in malignant strictures.

The Niti-S Fully Covered Biliary Stent is intended for maintaining biliary luminal patency in malignant and/or benign strictures.

The Bumpy Stent [TKxxyyCW or BKxxyyCW] is intended for use in malignant and/or benign biliary strictures and benign pancreatic ductal strictures.

WARRANTY

Taewoong Medical Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical

procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

4. Contraindication

The Niti-S & ComVi Biliary Stent is contraindicated for, but are not limited to:

- Patient with ascites.
- Placement in polypoid lesions.
- Patient with bleeding disorder.
- Intra-abdominal abscess.
- Patients with coagulopathy
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered/bare Stents is contraindicated. (see Warnings).
- Suspected or impending perforation.
- Patients for whom endoscopic techniques are contraindicated.
- Biliary obstruction preventing either endoscopic or percutaneous.
- Recapturing a stent during its deployment is contraindicated.

5. Warnings

- The safety and efficacy of this device for use in the vascular system has not been established.
- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration or/and fracture.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Niti-S Full Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- Niti-S Full Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Full Covered Stent can be repositioned immediately after deployment.
- Uncovered/bare Stents should not be removed once fully deployed; see Contraindications.
- Do not attempt to recapture/reload a stent once its deployment is advanced.

6. Potential complications

Potential complications associated with the use of Niti-S & ComVi Stent may include, but are not limited to:

Procedural Complications

- Bleeding
- Stent misplacement or inadequate expansion
- Pain
- Death (other than that due to normal disease progression)
- Intestinal perforation

Post Stent Placement Complications

- Bleedings
- Pain
- Perforation
- Bowel impaction
- Stent misplacement or migration
- Stent occlusion
- Tumor in-growth
- Tumor over-growth
- Stent fracture
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Acute cholecystitis
- Pancreatitis
- Cholangitis/Cholestasis
- Constipation
- Diarrhea
- Infection
- Liver abscess
- Peritonitis
- Sludge occlusion
- Ulcerations

7. Equipment required

● Percutaneous Placement

- 0.035” (0.89mm) guidewire at least 180cm long (preferably stiff or extra stiff)
- Introducer sheath appropriately sized for stent and introducer system (7 Fr or larger-uncovered 8Fr-covered)

● Endoscopic Placement

- 0.035” (0.89mm) guidewire at least 450cm long (preferably jag wire)
- Introducer sheath appropriately sized for stent and introducer system
- Endoscope system appropriately sized for instrument channel (8Fr or larger uncovered and covered, 3.7mm working channel)

● Short-wire introducer system

- Short-wire locking device
- Endoscope with Minimum 3.2mm working channel.
- 0.035”(0.89mm) stiff guidewire

8. Precautions

Read the entire User’s Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date “Use by”. Do not use the device beyond the use by date.
- The Niti-S & ComVi Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S & ComVi Stent is intended for single use only. Do not resterilize and/or reuse the device.

9. Instructions in the event of Damage

WARNING: Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

10. Procedure

- (50cm introducer usable length)

Percutaneous transhepatic cholangiography (PTC) should be performed prior to placement of the Niti-S & ComVi Stent to characterize the biliary tract morphology and extent of the malignant disease.

- (180cm introducer usable length)

Endoscopic retrograde cholangiopancreatography (ERCP) should be performed prior to placement of the Niti-S & ComVi Stent to characterize the biliary tract morphology and extent of the malignant disease.

① Examine stricture endoscopically and fluoroscopically

- a) Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.
- b) The Internal luminal diameter should be measured exactly with endoscope and/or fluoroscope.

② Stent Size Determination

- a) Measure the length of the target stricture.
- b) Select a stent size that is 20 to 40mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
- c) Measure the diameter of the reference stricture - it is necessary to select a stent which has an unconstrained diameter about 1 to 4mm larger than the largest reference target diameter, to achieve secure placement.

③ Stent Deployment Preparation

- The Niti-S & ComVi Stent can be placed with the aid of fluoroscopy and/or endoscopy.
- Pass a 0.035" (0.89 mm) guidewire to the level of the stricture.

A. Fluoroscopy Procedure

- a) Under the fluoroscopy guidance, insert a guide wire across the stricture to where the stent introducer system will be placed over the guide wire.
- b) Remove the stylet from the distal end of the introducer.
- c) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
- d) Flush the inner lumen of introducer system.

B. Endoscopy Procedure

- a) Under the endoscopic guidance, insert an endoscope to the level of the obstruction, then introduce the guide wire through the working channel of the endoscopy. Advance until the guide wire across the target stricture to where the stent introducer system will be placed over the guide wire.
- b) Remove the stylet from the distal end of the introducer.
- c) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
- d) Flush the inner lumen of introducer system.

C. Short-wire procedure

- a) Ensure that the Short-wire locking device (supplied separately) is firmly seated on the endoscope against the side of the working channel port as outlined in its direction for use.
- b) Place the 0.035 in (0.89mm) guidewire through the endoscope, into the ampulla and across the biliary stricture.

- c) A sphincterotomy and predilatation of the biliary stricture may be performed prior to stent implantation at the discretion of the physician.
- d) Remove safety cylinder and stylet before loading guidewire.
- e) Insert the trailing end of the guidewire through the tip of the delivery system and advance in short strokes until it exits the guidewire access port. The guidewire access port allows the guidewire to exit the delivery system approximately 30cm from the tip of the device. The guidewire then lies along the side of the delivery system.
- f) Advance the system over the guidewire. Unlock the wire to pass the stent into the scope.
- g) Lock the guidewire into place and continue advancing the system through the scope until the stent is in the desired location.

④ Stent Deployment Procedure


Figure 3

PRECAUTION: Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent.

- a) Under the fluoroscope and/or endoscopic guidance, position the introducer system to the center of the target stricture exactly.
- b) Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently slide the Y-connector back along the pusher towards the hub.
- d) When the center X-ray marker reaches the center of target stricture, continue pulling back on the Y-connector until the stent is fully deployed. (See figure 3, 4)


Figure 4

CAUTION Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage to bile duct.

The Y and T-Type stents may be deployed under fluoroscope or endoscopically, however, the T & Y-Type stents have larger interstices within the central mesh which are detailed by larger radiopaque markers. A second stent may be placed within these larger interstices in order to ensure coverage and drainage of two or more ducts. The central radiopaque marker should ideally be placed at the point of the bifurcation.

⑤ After Stent Deployment

- a) Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- b) Carefully remove the introducer system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent

- expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- c) Balloon dilatation inside the stent can be performed on demand.

11. Perform routine post implant procedures

- a) Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complications.

12. Instructions for removal of Niti-S Full Covered Stents (see Warnings)

Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the retrieval string and/or collapse the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

Caution: Do not allow excessive force to remove the stent as it may cause disconnect to the retrieval string.

To reposition a Niti-S & ComVi Full Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

Please note: the stent can only be repositioned and/or removed proximally.

Reuse Precaution Statement

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your Taewoong Medical Co., Ltd. representative. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and /or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage: Store at room temperature(10~40°C).

Disposal Requirements: The Niti-S & ComVi Stent's contained introducer system. At the end of its used device it must be properly disposed of in compliance with local or hospital regulations that it is appropriately packaged and secured.

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