CE
0120

Niti-S Nagi™ Stent

Authorized representative in Europe
Représentant autorisé en Europe
Bevollmächtigter in Europa
Rappresentante autorizzato per l'Europa
Representante autorizado en Europa
Avrupa yetkili temsilcisi
Autoriseret repræsentant i Europa
유럽 대리인

Consult instructions for use
Consulter les instructions d'utilisation
Gebrauchsanleitung beachten
Consultare le istruzioni per l'uso
Consulte las instrucciones de uso
Kullanım talimatlarına bakın
Konsultér brugsanvisning
사용자 지침 참조

Catalogue No.
N° de référence
Katalog-Nr.
N. di catalogo
N° de catálogo
Katalog No.
Referencenummer
형명

Attention, consult instructions for use
Attention, consulter les instructions d'utilisation
Achtung: Gebrauchsanleitung beachten
Attenzione, consultare le istruzioni per l'uso
Atención, consulte las instrucciones de uso
Dikkat! Kullanım talimatlarına bakın.
Vær opmærksom, konsultér brugsanvisning
주의, 사용상 지시 참조

Serial No.
N° de série
Serien-Nr.
N. Seriale
N° de serie
Seri No.
Seriennummer
일련번호

Date of Manufacture
Date de fabrication
Datum der Herstellung
Data di produzione
Fecha de fabricación
Üretim Tarihi
Fabrikationsdato
제조일자

Do not use if package is damaged
Ne pas utiliser si l'emballage est abîmé
Nicht verwenden, wenn die Verpackung beschädigt ist!
Non utilizzare la confezione se danneggiata
No usar si el paquete está dañado.
Ambalaj hasar görmüşse tekrar kullanmayın
Må ikke bruges, hvis emballagen er beskadiget
포장이 파손된 경우 사용금지

Temperature limitation
Limites de température
Temperaturvorgaben
Intervallo termico
Límites de temperatura
Isi sınırlaması
Temperaturgrænser
보관 온도

Manufacturer
Fabricant
Hersteller
Produttore
Fabricante
Üretici
Fabrikant
제조사

Sterilized using ethylene oxide
Stérilisé à l'aide d'oxyde d'éthylène
Mit Ethylenoxid sterilisiert
Sterilizzato con ossido di etilene
Esterilizado por óxido de etileno
Etilen oksit kullanılarak sterilize edilmiştir
Steriliseret med ethylenoxid
E.O. Gas 멸균

Use by (Expiration Date)
A utiliser avant (date d'expiration)
Verwendbar bis (Verfallsdatum)
Utilizzare entro (data di scadenza)
Utilizar antes de (fecha de vencimiento)
Son Kullanma Tarihi
Bruges før (udløbsdato)
유효기간

Do not reuse
Ne pas réutiliser
Nicht wiederverwenden
Non riutilizzare
No reutilizar
Yeniden kullanmayın
Må ikke genbruges
재사용 금지

Do not resterilize
Ne pas stériliser à nouveau
Nur einmal sterilisieren!
Non ristilizzare
No reesterilizar
Yeniden sterilize etmeyin
Må ikke reesteriliseres
재멸균 금지

MR Conditional
RM conditionnelle
Bedingte magnetische Resonanzsicherheit
RM compatible
RM condicional
Koşullu MR
MR-betinget
MR 조건부 인증

TaeWoong
MEDICAL

MRI Information

MR Conditional

Niti-S Nagi™ Stent was determined to be MR-conditional. Non-clinical testing demonstrated that the Niti-S Nagi™ Stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, Niti-S Nagi™ Stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position Nagi™ Stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm ²	676-mm ²	7,984-mm ²	977-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

User's Manual

1. Description

The Niti-S Nagi™ Stent consists of an implantable metallic stent and an introducer system.

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Niti-S Nagi™ Stent

Figure 1. Stent Model

The Stent is loaded in the introducer system and upon deployment the body of the stent imparts an outward radial force while both heads impart a horizontal inward force, allowing transmural drainage.

Niti-S Nagi™ Stent is Fully Silicone Covered and can be removed; (see Warnings).


Figure 2. Introducer System (Endoscopic)

- The endoscopic introducer system has a usable length of 180cm

Endoscopic Type is recommended

- When approached endoscopically

2. Principle of Operation

The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher towards the hub. Retraction of the outer sheath releases the stent.

3. Indication for Use

The Niti-S Nagi™ Stent is intended for drainage of a pancreatic pseudocyst through a transgastric or transduodenal approach.

WARRANTY

Taewoong Medical Co., Ltd. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

4. Contraindication

The Niti-S Nagi™ Stent is contraindicated for, but are not limited to:

- Hemodynamic instability
- Severe coagulopathy
- All others than indication for use
- Recapturing a stent during its deployment is contraindicated.

5. Warnings

- The safety and efficacy of this device for use in the vascular system has not been established.
- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Niti-S Nagi™ Stent is intended for short-term implantation, less than 30 days. Long-term patency of this stent has not been established and regular observation is recommended.
- The stent cannot be resheathed once deployment has been initiated.
- Do not attempt to recapture/reload a stent once its deployment is advanced.

6. Precautions

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The Niti-S Nagi™ Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S Nagi™ Stent is intended for single use only. Do not resterilize and/or reuse the device.

7. Instructions in the Event of Damage

WARNING: Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

8. Potential Complications

Potential complications associated with the use of Niti-S Nagi™ Stent may include, but are not limited to:

Procedural Complications

- Bleeding
- Stent misplacement
- Inadequate expansion
- Migration
- Pain
- Perforation

Post Stent Placement Complications

- Bleeding
- Pain
- Perforation
- Stent dislocation
- Stent misplacement or migration
- Stent occlusion
- Fever
- Pancreatitis
- Abscess formation
- Haemorrhage
- Vomiting
- Pneumoperitoneum
- Intraperitoneal leakage
- Peritonitis
- Hematoma
- Inflammation or Infection
- Fistula
- Ulceration
- Sepsis
- Rupture of intracystic artery
- Stent removal failure

9. Equipments Required

- Percutaneous Placement
Percutaneous placement has not been established yet.
- Endoscopic Placement
 - 0.035" (0.89mm) guidewire at least 450cm long
 - Echoendoscope or therapeutic duodenoscope with instrument channel of 3.7mm or larger
 - Needle
 - Needle knife and/or balloon catheter
 - Stent
 - Forceps and/or Snare

10. Procedure

- **(50cm introducer usable length)**
Percutaneous procedure has not been established yet.
- **(180cm introducer usable length)**
Endoscopy Ultrasound (EUS) should be performed prior to placement of the Niti-S Nagi™ Stent to characterize the extent of the lesion and its morphology.

A. Stent Size Determination

- a) The stent length and diameter should be determined by the doctor after endoscopic and/or fluoroscopic inspection of the lesion.
- b) The length and diameter of the stent should be chosen in a way that both walls of the created transmural drainage remain tightly together, thus preventing any migration.

B. Stent Deployment Preparation

A. Fluoroscopic Procedure

Fluoroscopic procedure has not been established yet.

B. Endoscopic Procedure

- a) Under endoscopic guidance, insert an endoscope until reaching the lesion. Then introduce a needle through the working channel of the endoscope. Advance the needle until reaching the lesion and puncture it.
- b) After puncture of the lesion, insert a guide wire through the needle and advance it across the lesion. Remove the needle slowly and carefully.
- c) After needle removal, insert a dilating device such as a needle knife and/or balloon catheter along the guide wire until reaching across the lesion and dilate it.

- d) After dilating, remove carefully a dilating device
- e) Remove the stylet from the distal end of the introducer.
- f) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
- g) Flush the inner lumen of introducer system.


Figure 3

C. Stent Deployment Procedure

PRECUATION: Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent.

- a) Under the fluoroscope and/or endoscopic guidance, position the introducer system
- b) Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently slide the Y-connector back along the pusher towards the hub
- d) Under EUS and/or fluoroscopic guidance, deploy the distal flare of the stent, and then deploy the proximal flare of the stent under endoscopic guidance while making sure that the stent connects both walls together.


Figure 4

CAUTION Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage.

D. After Stent Deployment

- a) Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- b) Carefully remove the introducer system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- c) Balloon dilatation inside the stent can be performed if judged necessary.

11. Perform Routine Post-implant Procedures

- a) Assess the good position of the stent and effective drainage. A Stent may require up to 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complications.

12. Instructions for Removal of Niti-S Nagi™ Stent

(see Warnings)

Grasp the retrieval string with forceps and/or collapse the proximal end of the Stent with snare then carefully and gently retrieve the Stent. Since long-term patency has not been established, we recommend to remove within 30 days after implantation.

Caution: Do not apply excessive force to remove the stent as it may disconnect the retrieval string.

To reposition a Niti-S Nagi™ Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

Reuse Precaution Statement

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your representative of Taewoong Medical Co., Ltd. For a single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage: Store at room temperature (10~40°C).

Disposal Requirements: The Niti-S Nagi™ Stent's contained introducer system. After use, it must be properly disposed of in compliance with local or hospital regulations that it is appropriately packaged and secured.

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