# Niti-S & ComVi Enteral Colonic Stent

**Authorized representative in Europe**  
Représentant autorisé en Europe  
Bevollmächtigter in Europa  
Rappresentante autorizzato per l'Europa  
Representante autorizado en Europa  
Avrupa yetkili temsilcisi  
Autoriseret repræsentant i Europa  
유럽 대리인

**Consult instructions for use**  
Consulter les instructions d'utilisation  
Gebrauchsanleitung beachten  
Consultare le istruzioni per l'uso  
Consulte las instrucciones de uso  
Kullanım talimatlarına bakın  
Konsultér brugsanvisning

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**Catalogue No.**  
N° de référence  
Katalog-Nr.  
N. di catalogo  
N° de catálogo  
Katalog No.  
Referenzenummer  
형명

**Attention, consult instructions for use**  
Attention, consulter les instructions d'utilisation  
Achtung: Gebrauchsanleitung beachten  
Attenzione, consultare le istruzioni per l'uso  
Atención, consulte las instrucciones de uso  
Dikkat! Kullanım talimatlarına bakın.  
Vær opmærksom, konsultér brugsanvisning  
주의, 사용상 지시 참조

**Serial No.**  
N° de série  
Serien-Nr.  
N. Seriale  
N° de serie  
Seri No.  
Serienummer  
일련번호

**Date of Manufacture**  
Date de fabrication  
Datum der Herstellung  
Data di produzione  
Fecha de fabricación  
Üretim Tarihi  
Fabrikationsdato  
제조일자

**Do not use if package is damaged**  
Ne pas utiliser si l'emballage est abîmé  
Nicht verwenden, wenn die Verpackung beschädigt ist!  
Non utilizzare la confezione se danneggiata  
No usar si el paquete está dañado.  
Ambalaj hasar görmüşse tekrar kullanmayın  
Må ikke bruges, hvis emballagen er beskadiget  
포장이 파손된 경우 사용금지

**Temperature limitation**  
Limites de température  
Temperaturvorgaben  
Intervallo termico  
Límites de temperatura  
Isi sınırlaması  
Temperaturgrænser  
보관 온도

**Manufacturer**  
Fabricant  
Hersteller  
Produttore  
Fabricante  
Üretici  
Fabrikant  
제조사

**Sterilized using ethylene oxide**  
Stérilisé à l'aide d'oxyde d'éthylène  
Mit Ethylenoxid sterilisiert  
Sterilizzato con ossido di etilene  
Esterilizado por óxido de etileno  
Etilen oksit kullanılarak sterilize edilmiştir  
Steriliseret med ethylenoxid  
E.O Gas 멸균

**Use by (Expiration Date)**  
A utiliser avant (date d'expiration)  
Verwendbar bis (Verfallsdatum)  
Utilizzare entro (data di scadenza)  
Utilizar antes de (fecha de vencimiento)  
Son Kullanma Tarihi  
Bruges før (udløbsdato)  
유효기간

**Do not reuse**  
Ne pas réutiliser  
Nicht wiederverwenden  
Non riutilizzare  
No reutilizar  
Yeniden kullanmayın  
Må ikke genbruges  
재사용 금지

**Do not resterilize**  
Ne pas stériliser à nouveau  
Nur einmal sterilisieren!  
Non ristilizzare  
No reesterilizar  
Yeniden sterilize etmeyin  
Må ikke resteriliseres  
재멸균 금지

**MR Conditional**  
RM conditionnelle  
Bedingte magnetische Resonanzsicherheit  
RM compatible  
RM condicional  
Kosullu MR  
MR-betinget  
MR 조건부 입증

**TaeWoong**  
MEDICAL

## ***MRI Information***

### **MR Conditional**

Niti-S & comvi stent was determined to be MR-conditional. Non-clinical testing demonstrated that the Niti-S & comvi stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

#### ***Static Magnetic Field***

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

#### ***MRI-Related Heating***

In non-clinical testing, Niti-S & comvi stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

#### ***Artifact Information***

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position Niti-S & comvi stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm <sup>2</sup>	676-mm <sup>2</sup>	7,984-mm <sup>2</sup>	977-mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

## User's Manual

### 1. Description

The Niti-S & ComVi Enteral Colonic Stent consists of the implantable metallic stent and introducer system.

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Niti-S Enteral Colonic Uncovered Stent
Niti-S Enteral Colonic Covered Stent
ComVi Enteral Colonic Stent


**Figure1. Stent Model**

The Stent is loaded into an introducer system and upon deployment the stent imparts an outward radial force on the luminal surface of the colorectum to establish patency. There are two types of introducer; TTS and OTW. (See figure 2, 3)

Niti-S Full Covered Enteral Colonic Stents used in benign stricture can be removed; (see Warnings).


Full Covered Enteral Colonic stents can be repositioned after deployment (see Warnings).

#### A. TTS Type Introducer

**Figure 2. TTS Type Introducer System**

#### B. OTW Type Introducer

**Figure 3. OTW Type Introducer System  
(Distal release & Proximal release)**

TTS means **T**hrough **T**he endoS**C**ope, and OTW means **O**ver **T**he **W**ire. The TTS type may be utilised for endoscopic delivery as the profile of the introducer system is compatible with a minimum 3.7mm working channel. The OTW type is not compatible with the endoscope due to the shorter length of the introducer system.

**TTS** type is recommended

- Where the target stricture region is narrow
- Where long-distance deploying is required
- Where exact and precise deployment is required

**OTW** is recommended

- Where  $\varnothing$  18, 20, 24, 26, 28mm stent is required
- Where fluoroscope is used

### 2. Principle of Operation

- For OTW Type; Distal release & TTS Type Introducer systems, The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher (2<sup>nd</sup> inner catheter) towards the hub. Retraction of the outer sheath releases the stent.

- For OTW Type; Proximal release systems, the distal part of outer sheath is pushed forward by immobilizing the Connector in one hand, grasping the hub with the other hand, and gently sliding the hub along the 2<sup>nd</sup> inner catheter towards the Connector. Retraction of the outer sheath releases the stent

### **3. Indication for Use**

The Niti-S & ComVi Enteral Colonic Stent is intended for maintaining colon luminal patency in colon stricture caused by intrinsic and/or extrinsic malignant and/or benign stricture.

### **WARRANTY**

Taewoong Medical Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the repair or replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

### **4. Contraindication**

The Niti-S & ComVi Enteral Colonic Stent is contraindicated for, but not limited to;

- Patient with ascites.
- Placement in polypoid lesions.
- Patient with bleeding disorder.
- Intra-abdominal abscess.
- Patients with coagulopathy.
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered/bare Stents is contraindicated. (See Warnings).
- Suspected or impending perforation.
- Enteral ischemia.
- Recapturing a stent during its deployment is contraindicated.

### **5. Warnings**

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Niti-S Full Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.

- Niti-S Full Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Full Covered Stent can be repositioned immediately after deployment.
- Uncovered/bare Stents should not be removed once fully deployed, (See Contraindications).
- Do not attempt to recapture/reload a stent once its deployment is advanced.

## **6. Potential complications**

Potential complications associated with the use of the colon stent may include, but are not limited to:

### Procedural Complications

- Bleeding
- Stent misplacement or inadequate expansion
- Pain
- Death (other than that due to normal disease progression)
- Intestinal Perforation

### Post Stent Placement Complications

- Bleeding
- Pain
- Perforation
- Bowel impaction
- Stent misplacement or migration
- Stent occlusion
- Tumor overgrowth
- Tumor ingrowth
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Constipation
- Diarrhea
- Infection
- Peritonitis
- Ulceration
- Symptoms of tenesmus or urgency/incontinence

## **7. Equipment required**

- Fluoroscope and/or Endoscope (with a channel size of 3.7 mm or greater)
- 0.035 in/ 0.89 mm guidewire (TTS)
- 0.038 in/ 0.97 mm guidewire (OTW; Distal & Proximal release)
- Introducer sheath appropriately sized for stent and introducer system

## **8. Precautions**

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.

- Check the expiration date “Use by”. Do not use the device beyond the use by date.
- The Niti-S & ComVi Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S & ComVi Stent is intended for single use only. Do not resterilize and/or reuse the device.

## **9. Instructions in the event of Damage**

**WARNING:** Visually inspect the system for any sign of damage. **DO NOT USE** if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

## **10. Procedure**

### **① Examine stricture endoscopically and fluoroscopically**

- a) Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.
- b) The Internal luminal diameter should be measured exactly with endoscope and/or fluoroscope.

### **② Stent Size Determination**

- a) Measure the length of the target stricture.
- b) Select a stent size that is 20 to 40mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
- c) Measure the diameter of the reference stricture - it is necessary to select a stent which has an unconstrained diameter about 1 to 4mm larger than the largest reference target diameter, to achieve secure placement.

### **③ Stent Deployment Preparation**

- The Niti-S & ComVi Stent can be placed with the aid of fluoroscopy, and/or endoscopy.
- Pass a 0.035” (0.89 mm) guidewire to the level of the stricture. (TTS)
- Pass a 0.038” (0.97 mm) guidewire to the level of the stricture. (OTW)

#### **A. TTS procedure**

- a) Under the endoscopic guidance, insert an endoscope to the level of the obstruction, then introduce the guide wire through the working channel of the endoscopy. Advance until the guide wire across the target stricture to where the stent introducer system will be placed over the guide wire.
- b) Remove the stylet from the distal end of the introducer.
- c) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
- d) Flush the inner lumen of introducer system.

#### **B. OTW procedure**

- a) Remove the stylet from the distal end of the introducer.
- b) Under the fluoroscopy guidance, insert a guide wire across the stricture to where the stent introducer system will be placed over the guide wire.
- c) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.

#### ④ Stent Deployment Procedure


Figure 4

**PRECAUTION:** Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent

#### A. OTW(Distal Release System) & TTS Introducer System

- Under the fluoroscope and/or endoscopic guidance, position the introducer system to the center of the target stricture exactly.
- Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in an anti-clockwise direction.
- To begin stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently slide the Y-connector back along the pusher towards the hub.
- When the center X-ray marker reaches the center of target stricture, continue pulling back on the Y-connector until the stent is fully deployed. (See figure 4, 5)


Figure 5

**CAUTION** Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage to the target or stricture.

#### B. OTW(Proximal Release System)

- Under the fluoroscope and/or endoscopy guidance, position the introducer system to the center of target stricture exactly.
- Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Connector by turning the valve more than twice in an anti-clockwise direction. The stent is now ready for deployment
- To begin stent deployment, immobilize the Connector in one hand and grasp the hub with the other hand. Gently slide the hub forward along the 2<sup>nd</sup> inner catheter towards the Connector.
- When the center X-ray marker reaches the center of target stricture, continue forwarding toward the connector until the stent is fully deployed. (See figure 4, 6)


Figure 6

**CAUTION** Do not push forward or pull backward on the Connector with the stent partially deployed. The Connector must be securely immobilized. Inadvertent movement of the Connector may cause misalignment of the stent and possible damage to the colorectum.

#### **⑤ After Stent Deployment**

- a) Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- b) Carefully remove the introducer system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- c) Balloon dilatation inside the stent can be performed if on demand.

#### **11. Perform routine post implant procedure**

- a) Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complications.

#### **12. Instructions for removal of Niti-S Full Covered Stents (see Warnings)**

Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the retrieval string and/or collapse the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

Caution: Do not allow excessive force to remove the stent as it may cause disconnect to the retrieval string.

To reposition a Niti-S & ComVi Full Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

Please note: the stent can only be repositioned and/or removed proximally.

#### **Reuse Precaution Statement**

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your Taewoong Medical Co., Ltd. representative. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and /or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Storage:** Store at room temperature(10~40°C).

**Disposal Requirements:** The Niti-S & ComVi Stent's contained introducer system. At the end of its used device it must be properly disposed of in compliance with local or hospital regulations that it is appropriately packaged and secured.

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