



UVENTATM Urethral Stent

Self-Expandable Metallic Stent



*"Minimal invasive treatment of urethral obstruction
for patient's quality of life."*

■ Fluoroscopic images



MANUFACTURER : TAEWOONGMEDICAL Co.,Ltd 14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do,
Korea 415-873 TEL:+82.31.996.0641~4 FAX:+82.31.996.0645~6

SALES OFFICE Bldg.C 4F, Daebangtriplaon BusinessTower, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do,
Korea 410-315 TEL:+82.31.904.6153 FAX:+82.31.904.6157

E-MAIL:sales@stent.net www.stent.net


 **TaeWoongMedical**

UVENTA™ Urethral Stent

● Features & Benefits

● Structure


- Different magnitudes of segmental radial force prevents complication such as reactive hyperplasia
- Full silicone covered to prevent tissue in-growth



● Radiopaque marker

3(three) at both ends

* CE is pending now



● Anti-migration flaps

4(four) flaps at the distal end prevent the risk of migration



● Removal string

Retrieval string at the proximal end facilitates easy removal of stent

● Ordering Information

Endoscopic Approach				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
US1004F	10	4	10	50
US1005F		5		
US1006F		6		
US1008F		8		
US1010F		10		
US1204F	12	4		
US1205F		5		
US1206F		6		
US1208F		8		
US1210F		10		
US1404F	14	4		
US1405F		5		
US1406F		6		
US1408F		8		
US1410F		10		
US1604F	16	4		
US1605F		5		
US1606F		6		
US1608F		8		
US1610F		10		

● Delivery system

