


# Bonopty<sup>®</sup>

Bone Biopsy System, 12 gauge




Optimize tactility and efficiency during  
your procedures with Bonopty<sup>®</sup> products

**AprioMed**

## Bonopty® Bone Biopsy System, 12 gauge

- New, larger gauge gains access into the bone, even through thick cortical bone
- Coaxial system facilitates multiple sampling, or treatment of lesions
- Achieve excellent core samples with few crush artifacts
- Successfully sample somewhat sclerotic lesions
- Secure sample in biopsy cannula during retrieval using Core Lock



Use of Bonopty® Penetration Set and Bonopty® Biopsy Set to sample a lytic lesion with bone structures in lumbar vertebrae.



Before using Bonopty® read the Instructions for Use which accompany the product. Bonopty® is a registered trademark of AprioMed AB, Uppsala, Sweden. Bonopty® is manufactured under one or more of the following patents: US 5,423,824; US 5,810,826; US 5,868,684; EP 0 794 733; JP 3710475; SE 470 177. Other patents pending.

### ORDERING INFORMATION

#### Bonopty® Bone Biopsy System, 12G

Product number	Product Description	Minimum Order Quantity
12-1272	Bonopty® Penetration Set, 12 gauge, 10 cm	5 units
12-1273	Bonopty® Biopsy Set, 13 gauge, 16.5 cm	5 units
12-1274	Bonopty® Extended Drill, 13 gauge, 16.5 cm	5 units

#### Bonopty® Bone Biopsy System, 14G

Product number	Product Description	Minimum Order Quantity
10-1062	Bonopty® Penetration Set, 14 gauge, 6.5 cm	5 units
10-1072	Bonopty® Penetration Set, 14 gauge, 9.5 cm	5 units
10-1082	Bonopty® Penetration Set, 14 gauge, 13 cm	5 units
10-1063	Bonopty® Biopsy Set, 15 gauge, 13 cm	5 units
10-1073	Bonopty® Biopsy Set, 15 gauge, 16 cm	5 units
10-1083	Bonopty® Biopsy Set, 15 gauge, 19.5 cm	5 units
10-1074	Bonopty® Extended Drill, 15 gauge, 16 cm	5 units

29157 Rev. 02 2012-12

Head Office:  
AprioMed AB  
Virdings Allé 28  
SE-754 50 Uppsala, Sweden  
Phone: +46 18 430 14 40  
Free Fax (Int'l): +46 770 339972

E-mail: [info@apriomed.com](mailto:info@apriomed.com)  
[www.apriomed.com](http://www.apriomed.com)

US Sales Office:  
AprioMed, Inc.  
2 Palmer Dr, Suite 1  
Londonderry, NH 03053, USA  
Phone: +1 603 421 0875  
Fax (US toll free): 866 239 1797