



# NAGI™ Stent

*for pancreatic pseudocyst drainage*



*"An endoscopic approach that uses this new FCSEMS (NAGI™) is feasible in the treatment of PFCs of both pancreatic pseudocysts and walled-off pancreatic necrosis."*

◀ Placement of a transnasal drainage tube for irrigation and direct endoscopic necrosectomy through the stent


by Naysuyo Yamamoto, MD [Gastrointest Endosc. 2013 May;77(5):809-14]

# Niti-S™ NAGI™ Stent


## for pancreatic pseudocyst drainage

### ● Features

- Wide and smooth flare edges
  - Prevent the risk of migration and possibility of stent related luminal damages
- Available in various diameters (Up to 16mm)
  - Optimize drainage and provide enough path for following necrosectomy in case
- Retrieval String for repositioning or easy removal
- Radiopaque marker : 3(three) at both ends & 2(two) in the middle



A large amount of pus came out through the stent



Conventional upper GI endoscope advances into the pseudocyst through the stent. Then performing of necrosectomy is using the snare forceps

[J Hepatobiliary Pancreat Sci. 2013 Mar;20(3):403-6]

### ● Ordering Information

Endoscopic Approach				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS1001FW	10	1	9	180
BS1002FW		2		
BS1003FW		3		
BS1201FW	12	1	10	
BS1202FW		2		
BS1203FW		3		
BS1402FW	14	2	10	
BS1403FW		3		
BS1602FW	16	2	10	
BS1603FW		3		

### Released Articles

- \* *Clinical evaluation of endoscopic ultrasonography-guided drainage using a novel flared-type biflanged metal stent for pancreatic fluid collection* by Shuntaro Mukai et al [Endosc Ultrasound. 2015 Apr-Jun;4(2):120-5]
- \* *EUS-guided drainage of hepatic abscess and infected biloma using short and long metal stents (with videos)* by Ryosuke Tonozuka, MD et al [Gastrointest Endosc. 2015;81(6):1463-9]
- \* *First report of endoscopic ultrasound-guided cholecystogastrostomy with a Nagi covered metal stent for palliation of jaundice in extrahepatic biliary obstruction* by Shuntaro Mukai et al [Endosc Ultrasound. 2015 Apr-Jun;4(2):120-5]
- \* *Preliminary report on a new, fully covered, metal stent designed for the treatment of pancreatic fluid collections* by Natsuyo Yamamoto, MD et al [Gastrointest Endosc. 2013 May;77(5):809-14]
- \* *Treatment of infected pancreatic pseudocysts using a novel, dedicated covered self-expandable metal stent (CSEMS) with an effective antimigration system* by I. Tarantino et al [Endoscopy. 2012;44 Suppl 2 UCTN:E147-8]
- \* *New fully-covered self-expandable metal stent for endoscopic ultrasonography-guided intervention in infectious walled-off pancreatic necrosis (with video)* by Takao Itoi et al [J Hepatobiliary Pancreat Sci. 2013 Mar;20(3):403-6]

