


Biliary Stent (Covered)




Niti-S™ Biliary Stent (Covered)

for benign and malignant biliary strictures

Features

- Fixed cell with braided construction
 - Flexible and resistant to fracture
- Atraumatic ends
 - Less hyperplasia at the edges
- Retrieval string facilitates safe and smooth removal
- Silicone covering on both inner and outer surface
 - Prevent the risk of tumor ingrowth
 - Help smooth bile flow



- Radiopaque marker: 3 (three) at both ends & 2 (two) in the middle

Ordering Information

Endoscopic Approach					Percutaneous Approach				
► Fully covered					► Fully covered				
Code	Stent		Delivery		Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
BS0604F	6	4	8.5	180	TS0604F	6	4	8.5	50
BS0605F		5			TS0605F		5		
BS0606F		6			TS0606F		6		
BS0607F		7			TS0607F		7		
BS0608F		8			TS0608F		8		
BS0609F		9			TS0609F		9		
BS0610F		10			TS0610F		10		
BS0612F		12			TS0612F		12		
BS0804F		8			4		TS0804F		
BS0805F	5				TS0805F	5			
BS0806F	6				TS0806F	6			
BS0807F	7				TS0807F	7			
BS0808F	8		TS0808F	8					
BS0809F	9		TS0809F	9					
BS0810F	10	TS0810F	10						
BS0812F	12	TS0812F	12						
BS1004F	10	4	TS1004F	4					
BS1005F		5	TS1005F	5					
BS1006F		6	TS1006F	6					
BS1007F		7	TS1007F	7					
BS1008F		8	TS1008F	8					
BS1009F		9	TS1009F	9					
BS1010F		10	TS1010F	10					
BS1012F		12	TS1012F	12					

- * Available both ends 5mm bare [Coding: BS**__B / TS**__B]
- * Short-wire delivery system is available [Coding: BSM**__F / BSM**__B]

Released Article

- * Clinical Outcome of Endoscopic Ultrasound-Guided Liver Abscess Drainage Using Self-Expandable Covered Metallic Stent (with Video) by Takeshi Ogura et al [Dig Dis Sci. 2016 Jan;61(1):303-8]
- * Polyurethane-Covered Self-Expandable Nitinol Stent for Malignant Biliary Obstruction: Preliminary Results by Young-Min Han et al [Cardiovasc Intervent Radiol. 2002 Sep-Oct;25(5):381-7]