





DOUBLE™

Esophageal stent (Anti-reflux)



MANUFACTURER : TAEWOONG MEDICAL Co., Ltd 14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, South Korea 415-873 TEL: +82 31 996 0641~4 FAX: +82 31 996 0645~6

SALES OFFICE Bldg.C 4F, Daebangtriplaon BusinessTower, 158 Haneulmaeul-ro, ilsandong-gu, Goyang-si, Gyeonggi-do, South Korea 410-315 TEL: +82 31 904 6153 FAX: +82 31 904 6157

E-MAIL : sales@stent.net www.stent.net



Taewoong Niti-S™ Double™ Esophageal stent (Anti reflux) for preventing gastroesophageal reflux

● Features

- **Double layered design with Anti-reflux skirt**
 - PTFE skirt blocks gastric reflux with the stent placement at EG junction
 - Additional uncovered outer mesh helps to resist migration
- **Retrieval string** at proximal end helps repositioning
- **Radiopaque marker** : 4(four) at both ends & 2(two) in the middle

● Ordering Information

Distal Release					
Code	Stent		Delivery		
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)	
EA1606FD	16	6	18	70	
EA1608FD		8			
EA1610FD		10			
EA1612FD		12			
EA1614FD		14			
EA1615FD		15			
EA1806FD	18	6			
EA1808FD		8			
EA1810FD		10			
EA1812FD		12			
EA1814FD		14			
EA1815FD		15			
EA2006FD	20	6	20		
EA2008FD		8			
EA2010FD		10			
EA2012FD		12			
EA2014FD		14			
EA2015FD		15			
EA2206FD	22	6			22
EA2208FD		8			
EA2210FD		10			
EA2212FD		12			
EA2214FD		14			
EA2215FD		15			
EA2406FD	24	6	22		
EA2408FD		8			
EA2410FD		10			
EA2412FD		12			
EA2414FD		14			
EA2415FD		15			
EA2806FD	28	6			
EA2808FD		8			
EA2810FD		10			
EA2812FD		12			
EA2814FD		14			
EA2815FD		15			

